



Handicare 1000

A Minimum swivel radius – with arm spacers – with extended arm spacers	670mm 693mm 710mm
B Minimum wall to seat back	75mm†
C Seat depth – with extended seat	360mm 425mm
D Minimum folded width to footplate	385mm†
E Minimum open width to edge of footplate	650mm†
F Armrest width – external – with arm spacers – with extended arm spacers	550mm 590mm 637mm
G Armrest width – internal – with arm spacers – with extended arm spacers	460mm 516mm 562mm
H Seat back height	453mm
I Footplate to seat height – standard seat option – perch seat option^	465-540mm 677-877mm
J Minimum footplate height	85mm
Bottom parking space required	550mm (56°) min to 675mm at shallow angle (27°)††
Minimum track intrusion into staircase	115mm†
Minimum track overhang at top of stairs	86mm
Minimum track intrusion at bottom of stairs	400mm
Minimum track length with hinge	2300mm
Maximum weight limit – standard version – heavy duty version (XL)* – extra heavy duty version (XXL)*	137kg/21.5 stone 160kg/25 stone 200kg/31.5 stone
Speed	0.12m/sec (approx)
Power supply	24v DC battery
Charger	100-240v (50-60Hz)
Maximum staircase angle	56° (52° with hinge)

† +45mm for extra heavy duty (XXL) 2 way swivel with extended arm spacers.

†† +20mm with arm spacers, +44mm with extended arm spacers.

^Not available for extra heavy duty (XXL). *Not available with hinge.



Customer enquiries: 01823 661037
Freephone: 0800 6696 720

Please be aware that actual colours may vary slightly from those shown. All information is provided in good faith but does not constitute a guarantee. Handicare reserves the right to alter designs and specifications without prior notice.



handicare



Handicare 1000
Straight stairlift



The benefits of a Handicare 1000 straight stairlift

If you are thinking about installing a stairlift the Handicare 1000 offers you a safe and reliable solution.

Handicare stairlifts have helped thousands of people overcome the problems posed by going up and down their stairs at home.

The Handicare 1000 stairlift offers you all the basic functions that you would expect plus the option to upgrade to a powered seat swivel and/or folding footplate.

With a choice of six upholstery colours the Handicare 1000 will compliment most colour schemes.



The slimline track is fitted to the stairs, not the wall, so it is unlikely you will need to redecorate following the installation of your Handicare 1000 stairlift.

Like all Handicare stairlifts the 1000 Series has a battery backup which means it will continue to work in the event of a power cut so you do not need to worry about being stranded up or downstairs.

The compact unit folds up leaving the staircase free for other users and the slim track maximises clear space on the stairway and can be supplied with a manual or powered hinged track if required.

The Handicare 1000 also comes with two remote controls which means, once folded, you can easily send it to the upstairs charge point freeing up space in your hallway.



Safety as standard

Your Handicare 1000 stairlift will come with a range of safety sensors that are designed to stop your stairlift if it meets any obstruction as you travel up the stairs.

Your Handicare 1000 stairlift has an easy release retractable seat belt for added peace of mind. If you prefer you can have an easy fasten Velcro seat belt.

Easy to use

The operating toggle is easy to use, simply hold it towards the direction you wish to travel, up or down, and your Handicare 1000 will move smoothly and gently on its way.

The control is conveniently located at the end of the arm leaving both arms clear to provide unhindered support when getting on or off the seat.