



Handicare Simplicity+

A Minimum swivel radius – with arm spacers – with extended arm spacers	655mm 678mm 695mm
B Minimum wall to seat back	65mm
C Seat depth – with extended seat	360mm 425mm
D Minimum folded width to footplate	370mm
E Minimum open width to edge of footplate – large footplate (perch seat)	570mm 630mm
F Armrest width – external – with arm spacers – with extended arm spacers	550mm 590mm 637mm
G Armrest width – internal – with arm spacers – with extended arm spacers	460mm 516mm 562mm
H Seat back height	453mm
I Footplate to seat height – standard seat option – perch seat option	430-540mm 677-877mm
J Minimum footplate height	70mm
Bottom parking space required	575mm (50°) to 690mm at shallow angle (30°) [†]
Minimum track intrusion into staircase	185mm
Minimum track overhang at top of stairs (typical at 45°)	150mm
Minimum track intrusion at bottom of stairs	250mm
Maximum weight limit – Hinge/Slide Track	140kg/22 stone 120kg/19 stone
Speed	0.12m/sec (approx)
Power supply	24v DC battery
Charger	100-240v (50-60Hz)
Staircase angle*	27-50°

*Recommended minimum Slide Track staircase angle 37°.
[†] +20mm with arm spacers, +44mm with extended arm spacers.



Customer Enquiries: 01823 661037

Please be aware that actual colours may vary slightly from those shown. All information is provided in good faith but does not constitute a guarantee. Handicare reserves the right to alter designs and specifications without prior notice.



Handicare Simplicity+
 Straight stairlift



The benefits of a Handicare Simplicity+ straight stairlift

If you are thinking about installing a stairlift the Handicare Simplicity+ offers you a safe and reliable solution.

Handicare stairlifts have helped thousands of people overcome the problems posed by going up and down their stairs at home.

The Handicare Simplicity+ straight stairlift offers you all the basic functions that you would expect plus the option to upgrade to a powered seat swivel and/or folding footplate.

If you have an obstruction at the base of the stairs, you can choose from a manually operated hinged track or the unique Slide Track that is only available from Handicare.

With a choice of six upholstery colours the Handicare Simplicity+ will compliment most colour schemes.

The track is fitted to the stairs, not the wall, so it is unlikely you will need to redecorate following the installation of your Handicare Simplicity+ stairlift.

Like all Handicare stairlifts the Simplicity+ has a battery backup which means it will continue to work in the event of a power cut so you do not need to worry about being stranded up or downstairs.

The compact unit folds up leaving the staircase free for other users and the Handicare Simplicity+ also comes with two remote controls which means, once folded, you can move it along the track freeing up space in your hallway.



Safety as standard

Your Handicare Simplicity+ stairlift will come with a range of safety sensors that are designed to stop your stairlift if it meets any obstruction as you travel up the stairs.

Your Handicare Simplicity+ stairlift has an easy release retractable seat belt for added security and peace of mind. If you prefer you can have an easy fasten Velcro seat belt.

Easy to use

The operating toggle is easy to use, simply hold it towards the direction you wish to travel, up or down, and your Handicare Simplicity+ will move smoothly and gently on its way.

The control is conveniently located at the end of the arm leaving both arms clear to provide unhindered support when getting on or off the seat.

